

PRODUCT SPECIFICATION

Spec code: NOF-C0139
Product name: Nitrile Powder free
Date of issue: 17/06/14

GENERAL INFORMATION

Type: single-use disposable examination or protective glove
Marking: all information on dispenser box
Shape: ambidextrous – straight fingers
Material: Nitrile Butadiene Rubber (NBR)
Colour: white
Inner: online chlorinated
Outer: no treatment
Cuff/Surface: rolled rim/Fingertip textured

Thickness: <small>(double measured)</small>	Finger: min. 0.14 mm	Palm: min. 0.10 mm	Cuff: min. 0.08 mm
---	-----------------------------	---------------------------	---------------------------

Residual powder: <small>[according to EN 455-3]</small>	≤ 2 mg / glove	Force at break: <small>[according to EN 455-2]</small>	med. 6.0 N (during shelf life)
---	----------------	--	-----------------------------------

Length:	median 265 mm
----------------	---------------

Available Sizes:	XS (5-6)	S (6-7)	M (7-8)	L (8-9)	XL (9-10)
-------------------------	----------	---------	---------	---------	-----------

Storage: 3 years storage life (date of expiration see on the box). Keep area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolour the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

Quality certification: ISO 9001, ISO 13485
Conformity to guidelines: Personal Protective Equipment 89/686/EEC: CE cat. III
Medical Device Directive 93/42/EEC: CE Class I

Conformity to standards:

EN 420	EN 374-2: specified	EN 374-3: Class 6
EN 455	level 2 - AQL 1.5	Sodium Hydroxide 40%
ASTM D6319	level 3: AQL <0.65 (G1) level 2: AQL <1.5 (G1) level 1: AQL <4.0 (S4)	



Conformity to food regulation: 1935/2004 

WARNING AND REPORTING SYSTEM

Potentially allergenic ingredients: Dithiocarbamate type, Zink-2-Mercaptobenzothiazolate
(For further information, the actual version of the ingredients list is given on request.)

Medical Device Vigilance and Reporting System: According to the official reporting criteria of MDD, incidents caused by examination gloves or surgical gloves must be reported immediately to our medical device reporting officer. Tel.: +43 2630 310 ext. -510; Fax: -549

Caution! This product contains accelerators which may cause allergic reactions!

Established (PM/DM):	Reviewed (QA/DM):	Approved (HoM/DM):	Customer approval
 Mag.(FH) K. Miller	 DI Dr. René Stempffer	 Mag. B. Schnetzlinger	